

Using system functions

The System Functions window provides access to administrative tools and settings that control basic operating characteristics of the system. Among these features and options are:

- Shutting down/restarting the application software and the Max PC.
- Setting the unit's speaker volume.
- Enabling/disabling "on-demand" processing of Rx orders.
- Enabling/disabling Auto-Run and Auto-Retry operating modes.
- Enabling/disabling Auto-Run operating mode .
- Enabling/disabling the left and right vial feeders.
- Backing up and restoring cell configurations.
- Resetting the Max's robotic arm and running diagnostics.
- Using the Label Wizard.

System Functions window



The System Functions window is available on both sides of the Max. The version of this window displayed on the *Inventory side*, however, does not offer all the functions and options available on the *Prescription side*.

On the *Prescription side*, the System Functions window consists of 4 tabs.

The screenshot displays the Parata control panel interface. At the top left is the Parata logo with a lock icon. On the top right are icons for settings (gears) and home (house). Below these is a navigation bar with four tabs: General, Tools, Configuration, and Advanced. The main area is divided into several sections:

- Machine Functions:** Contains buttons for Restart, Close, and Shut Down.
- Volume Selection:** Features a volume slider with left and right arrow buttons, and labels for MUTE and MAX.
- Reports:** Includes a text description: "View reports on the machine's drug inventory and prescription history." and a Launch button.
- Diagnostic Functions:** Contains buttons for Home, Park, and Save Trace.
- Drug Finder:** Includes a text description: "Scan a stock bottle to view its drug parameters." and a Launch button.
- Scanner Echo:** Features a Format: label.


A large text box in the center contains the following information:

This machine is an automated system for labeling, filling, and capping prescription drugs. It is designed to process prescription orders rapidly with consistent counting accuracy. It easily adapts to most pill shapes and sizes and enables the customer to customize the system to accommodate individual pharmacy preferences and requirements.

Version: 4.5.0
Serial Number: 04B2808-10009

© 2007-2015 Parata Systems, LLC. All rights reserved.

At the bottom of the interface, there is a status bar showing a mail icon with the number 4, the text "Running", a dropdown arrow, and buttons for Run and Pause.



On the *Inventory side*, only the General tab functions (labeled General Functions) are available.

The image shows a control panel interface for a Parata machine. At the top left is the Parata logo with a lock icon. To the right are icons for settings (gears) and home (house). Below this is a 'General Functions' section with several rows of controls:

- Machine Functions:** Restart, Close, Shut Down
- Volume Selection:** A slider with left and right arrow buttons, labeled 'MUTE' and 'MAX'.
- Reports:** View reports on the machine's drug inventory and prescription history. Launch
- Diagnostic Functions:** Home, Park, Save Trace
- Drug Finder:** Scan a stock bottle to view its drug parameters. Launch
- Scanner Echo:** Format:


In the center, there is a text box with the following content:

This machine is an automated system for labeling, filling, and capping prescription drugs. It is designed to process prescription orders rapidly with consistent counting accuracy. It easily adapts to most pill shapes and sizes and enables the customer to customize the system to accommodate individual pharmacy preferences and requirements.

Version: 4.5.0
Serial Number: 04B2808-10009

© 2007-2015 Parata Systems, LLC. All rights reserved.

At the bottom left, there is a status indicator showing a box with the number '4' and the word 'Running'. To the right are buttons for Run and Pause.



Touch the System Functions button on either side of the machine to display the System Functions window.




General

The screenshot displays the Parata software interface. At the top left is the Parata logo with a lock icon. To the right are icons for settings and home. Below this is a navigation bar with four tabs: **General**, **Tools**, **Configuration**, and **Advanced**. The **General** tab is active and contains several sections:

- Machine Functions:** Includes buttons for **Restart**, **Close**, and **Shut Down**.
- Volume Selection:** Features a volume slider with **MUTE** and **MAX** labels, and left/right arrow buttons.
- Reports:** Includes a text description: "View reports on the machine's drug inventory and prescription history." and a **Launch** button.
- Diagnostic Functions:** Includes buttons for **Home**, **Park**, and **Save Trace**.
- Drug Finder:** Includes a text description: "Scan a stock bottle to view its drug parameters." and a **Launch** button.
- Scanner Echo:** Includes a **Format:** label.

A large text box in the center provides a description of the machine: "This machine is an automated system for labeling, filling, and capping prescription drugs. It is designed to process prescription orders rapidly with consistent counting accuracy. It easily adapts to most pill shapes and sizes and enables the customer to customize the system to accommodate individual pharmacy preferences and requirements." Below this text, it lists "Version: 4.5.0" and "Serial Number: 04B2808-10009", followed by a **Details** button. At the bottom of this box is the copyright notice: "© 2007-2015 Parata Systems, LLC. All rights reserved."

At the bottom of the interface is a status bar. On the left, it shows a status icon with the number "4" and the text "Running". On the right, there are buttons for **Run** and **Pause**, along with a dropdown menu icon.



From the General tab, you can perform the following tasks:


- Shut down/restart the application
- Shut down the Max.
- Adjust the speaker volume
- Create, display, export and print reports

- Save a diagnostic “trace” file.
- Reset mechanical subsystems
- Scan a stock bottle and view its drug parameters
- Create stock bottle mappings⁴
- Perform a “scanner echo”
- Check the software version and unit serial number

NOTE: On the *Inventory screen*, only the **General** tab is available. All other system functions must be performed on the *Prescription screen*. These functions, however, require the Parata Support password.

Shut down the system

From the **General** tab, you can shut down the Max applications, reboot the operating system and shut down the entire unit. These functions can be performed on either side of the Max.



- The **Restart** button closes and then restarts the Max applications.
- The **Close** button closes the Max applications only.
- The **Shut down** button shuts down the entire unit, including the PC.

The section [Error! Reference source not found.](#) provides more detailed information on shutdown procedures.

Diagnostic Functions: Home System



Touching this button returns all mechanical subsystems to their starting (“home”) positions.

NOTE: The **Home** button is accessible only when the Max is in Pause mode; otherwise, it is disabled.

⁴ Not shown. This option is displayed at sites using the DIN numbering method.

Diagnostic Functions: Park System



Touching this button moves all mechanical subsystems to safe positions in preparation for system maintenance.

NOTE: The Park button is accessible only when the Max is in Pause mode; otherwise, it is disabled.

Figure 25: The Home and Park buttons are disabled when the Max is in Run mode.



Diagnostic Functions: Save Trace



A "trace file" is a log of events occurring on the Max, used primarily by system technicians for diagnostic purposes.

Drug Finder

This utility lets you scan a stock bottle and view all information the MDL contains about that drug. If the drug is not in the MDL, a message will appear.



NOTE: The drug does not need to be housed in the Max.

When you touch the Launch button, you will be prompted to scan the bottle.



Either the information will display, or you will see a message that the drug is not recognized.

Figure 26: Drug information


Drug Parameters	
Drug NDC/DIN	59762-5012-01
Brand Name	Spironolactone 50 mg Tablet
Generic Name	Spironolactone 50 mg Tablet
Custom Name	
30 Dram Capacity	105
Nozzle Height	3.50
Nozzle Width	B.75
Baffle	4
Pressure	40 psi
Manufacturer	GREENSTONE LTD
Schedule	RX
Front Imprint	G
Back Imprint	5012
Cleaning Code	B
Recommended	Yes
Pre-pack Quantities	
Drug Code Mappings	

Figure 27: Drug not found

Drug Parameters	
Drug NDC/DIN	
Brand Name	Unrecognized Drug
Generic Name	Spironolactone 50 mg Table
Custom Name	
30 Dram Capacity	
Nozzle Height	
Nozzle Width	
Baffle	
Pressure	
Manufacturer	
Schedule	
Front Imprint	
Back Imprint	
Cleaning Code	
Recommended	
Pre-pack Quantities	
Drug Code Mappings	


Stock Bottle Mappings

This button, available *only* at installations that use the DIN numbering system, launches the Stock Bottle Mappings Wizard.



Scanner Echo


When the General tab is selected and you use either scanner to scan a bar code, the bar code's format and encoded value are displayed in the Scanner Echo field, as shown in the following examples:



About

The About dialog box (appearing below the Scanner Echo field) displays descriptive information about the Max, as well as the unit’s serial number and the version of the application software currently installed.

Figure 28: This window is similar to the one displayed on your system.



Touch the Details button to display additional unit-specific information.




Configuration


(Prescription side only)

NOTE: Access to these functions requires the Parata Support password. Please contact Parata Technical Assistance Center if you need a configuration item changed.

The options available on the Configuration tab let you customize certain elements of Rx order processing.



For example, Fill To Empty mode is enabled unless you elect to disable it.



From the Configuration tab, you can:

- Enable/disable Auto Run mode (default = disabled)
- Enable/disable Auto Rerun mode (default = *enabled*)
- Enable/disable Fill On Demand mode (default = disabled)
- Enable/disable Fill To Empty mode (default = *enabled*)
- Enforce/ignore package size (default = disabled)
- Enforce/ignore pre-pack quantities (default = disabled)
- Enable/disable Right Vial Bin use (default = *enabled*)
- Enable/disable Left Vial Bin use (default = *enabled*)
- Enable/disable Vial Labeling mode (default = *enabled*)
- Enable/disable Back-up Vial mode (default = *enabled*)


- Enable/disable Complete On Ready mode (default = disabled)

The Configuration tab displays “toggle” buttons that display the option’s current setting and let you change it.

A green checkmark indicates the option is enabled. When you touch this button, the option is *disabled*.

A red “X” indicates the option is disabled. When you touch this button, the option is *enabled*.


Label Wizard




The Label Wizard is only accessible to administrative / support users. Please contact Parata Technical Assistance Center if you need your label edited.

The Label Wizard is used to create and print cell labels and custom vial labels.

Figure 30: Label Wizard Main window






Advanced

(Prescription side only)

NOTE: In general, only Parata technicians should use these functions. With the exception of Audit Configuration, access to Advanced functions requires the Parata Support password.

From the Advanced tab you can access the following features and options:



- Advanced Diagnostics
- Label Wizard
- Prescription Routing
- Back Up Database
- Clean Database
- Machine Configuration
- Audit Configuration
- Fleet Management Configuration
- Camera Diagnostics

Backups

Your system's database is backed up automatically at predesignated intervals. These backups do not require operator participation—they run undetected.

NOTE: By default, script data retention for return to stocks older than 30 days is deleted.